

# INOGENONE® G3

## Portable Oxygen Concentrator

The INOGENONE® G3 is a popular choice amongst oxygen users for its ease of use, stylish and non-medical feel.

With up to 5 pulse flow settings available, the INOGENONE® G3 strikes the perfect balance between those wanting a lightweight solution and the flexibility of delivering oxygen at a more demanding level.



**Pulse Flow**

**1 to 5**  
Pulse flow settings



**Weight**

From  
**4.9lbs**  
(2.2kgs)



**Battery life**

Up to  
**8**  
Hours  
on setting 2\*



*intermedical*  
H O M E C A R E

[www.portableoxygen.co.uk](http://www.portableoxygen.co.uk) | 0800 028 2194

Trusted in Respiratory Health

# Smart, stylish design with the user in mind.

This compact oxygen solution weighs from only 4.9lbs (2.2kgs) and can be worn easily over the shoulder in its protective black bag. A backpack accessory is also available separately for even greater comfort and independence.

The G3 is powered by either an 8 or 16 cell battery that slots neatly onto the base of the machine with an approximate run-time of 4 or 8 hours respectively (on setting 2).

The user-interface is clear and easy to operate with a white back-lit LCD display that shows the remaining battery operation time in hours/minutes and the current setting.




◀ Backpack accessory is also available for the G3.



◀ The INOGENONE® G3 has two battery options: 8 cell & 16 cell. The battery re-charges when installed and the machine is connected to mains or DC power. A desktop charger is also available to charge any supplementary batteries. (optional accessory)

## 8 Cell Battery Package

- INOGENONE® G3 main unit
- 8 cell rechargeable battery
- Carry bag with shoulder straps
- AC mains power supply unit
- DC power supply unit
- User manual
- Quick start guide
- Nasal cannula



21.0 (h) x 22.2 (w) x 7.6 (d)  
in cm with 8 cell battery installed

## 16 Cell Battery Package

As above including:

- 16 cell rechargeable battery instead of 8 cell battery

### Approximate battery run-time (in hours:minutes)

	Standard Battery	Double Battery
Setting 1	4:30	9:00
Setting 2	4:00	8:00
Setting 3	3:00	6:00
Setting 4	2:00	4:00
Setting 5	1:40	3:20

\*Run-time dependent on breathing rate

Concentration	90% -3% /+ 6%
Weight	4.9lbs (2.2kg) - 8 cell batt. 5.7lbs (2.6kg) - 16 cell batt.

Power AC: 100-240v DC: 13.5-15.5v

Noise 39db(A) on setting 2

Batt. charge time 8 cell: 4 hours  
16 cell: 8 hours

Warranty 3 years on unit  
1 year on batteries, sieve beds & accessories

✈ Airline approved Yes, subject to airline's terms & conditions

~~CPAP~~  
~~BiPAP~~ This device may be suitable for use whilst sleeping. This device cannot be used with a CPAP or BiPAP device.