

# CAIRE® Freestyle Comfort

## Portable Oxygen Concentrator

The FreeStyle Comfort by CAIRE (Formerly AirSep) is a lightweight, multi-setting portable oxygen concentrator. Weighing from as little as 5lbs (2.3kgs), this stylish and ergonomically designed concentrator provides up to 5 pulse-flow settings.

Its unique curved design is aimed to conform to the user's body to reduce bouncing whilst walking and increase overall comfort when using the concentrator.



Pulse Flow

**1 to 5**  
Pulse flow  
settings



Weight

From  
**5.0lbs**  
(2.3kgs)



Battery life

Up to  
**8:30**  
Hours  
on setting 2\*



*intermedical*  
H O M E C A R E

[www.portableoxygen.co.uk](http://www.portableoxygen.co.uk) | 0800 028 2194

Trusted in Respiratory Health

# Comfort by name. Comfort by nature.

CAIRE®  
Freestyle Comfort

The FreeStyle Comfort is incredibly simple and easy to use with just three buttons to power on/off the machine and adjust the flow to the desired setting - The only buttons an oxygen user is ever likely to need on a portable oxygen concentrator.

The large LCD screen shows current setting and remaining battery life in colour. Should the Comfort sound an alarm i.e. low battery or detached nasal cannula, the user will be prompted with an easy to understand graphic on the display to help troubleshoot the problem.

The FreeStyle Comfort is available with a detachable 8 or 16 cell battery that conveniently slots onto the bottom of the machine.

It comes complete with a generous 3-year warranty on the main unit and sieve beds. Batteries and other accessories come with a standard 1-year warranty.




Two battery options: Get up to 8 hours run-time with an 8 cell battery or up to 16 hours run-time with a 16 cell battery.

The colour LCD display shows current setting and remaining battery life in %



## 8 Cell Battery Package

- CAIRE® FreeStyle Comfort main unit
- 8 cell rechargeable battery
- Carry bag
- Shoulder straps
- AC mains power supply unit
- DC power cable
- User manual
- Nasal cannula



## 16 Cell Battery Package

As above including:

- 16 cell rechargeable battery instead of 8 cell battery

25.4 (h) x 18.5 (w) x 7.9 (d)  
in cm with 8 cell battery inserted

27.9 (h) x 18.5 (w) x 7.9 (d)  
in cm with 16 cell battery inserted

Approximate battery run-time  
(in hours:minutes)

	8 Cell Batt.	16 Cell Batt.
Setting 1	8:00	16:00
Setting 2	4:15	8:30
Setting 3	3:00	6:00
Setting 4	2:15	4:30
Setting 5	2:00	4:00

\*Run-time dependent on breathing rate

Concentration	90% (-3% / + 5.5%)
Weight	5.0lbs (2.3kg) w/ 8 cell battery 6.0lbs (2.7kg) w/16 cell battery

Power	AC: 100-240v DC: 11-18v
Noise	43db(A) on setting 2
Batt. charge time	8 cell: Approx. 3.5 hours 16 cell: Approx. 6 hours
Warranty	3 years on unit & sieve beds 1 year on batteries & accessories
Airline approved	Yes, subject to airline's terms & conditions



This device may be suitable for use whilst sleeping. This device cannot be used with a CPAP or BiPAP device.

intermedical  
HOMECARE

0800 028 2194  
www.portableoxygen.co.uk

Intermedical (UK) Limited, Homecare Division  
Unit 6 Mill Hall Business Estate, Aylesford, Kent, ME20 7JZ, United Kingdom

CAIRE®  
UK AUTHORISED  
DISTRIBUTOR

globalgroup  
ISO 13485  
UKAS  
MANAGEMENT  
SYSTEMS  
5965  
Cert: 87M18032

Trusted in **Respiratory Health**