

PHILIPS

RESPIRONICS

SimplyGo Mini

Portable Oxygen Concentrator

The SimplyGo Mini is made by none other than one of the world's leading and trusted consumer brand. Philips Respironics have spared no expense in creating a portable oxygen concentrator that combines independence with contemporary living, with a completely non-medical appearance and feel.

Weighing as little as 5lbs, the SimplyGo Mini provides up to 5 pulse flow oxygen settings. This small and lightweight oxygen concentrator is easy to carry and comes in a stylish brown carry bag with shoulder strap. For additional comfort, the SimplyGo Mini can also be worn over both shoulders with the stylish backpack option (available separately).



Pulse Flow

1 to 5
Pulse flow
settings



Weight

From
5.0lbs
(2.3kgs)



Battery life

Up to
9
Hours
on setting 2*



intermedical
H O M E C A R E

www.portableoxygen.co.uk | 0800 028 2194

Trusted in Respiratory Health

Powerful and stylish. With touchscreen control.

Whilst mobile, the SimplyGo Mini is powered with a standard or extended rechargeable battery that slots neatly on to the base of the machine - with up to 4½ or 9 hours operation time respectively (based on setting 2). A car charger is also provided keeping users mobile for longer.

The user-friendly colour touchscreen allows for easy operation from the fingertip as well as displaying current setting and remaining battery life. SimplyGo Mini is user maintenance-free with no filters to clean or replace. Accessory bag provided as standard to keep all power cables together.



◀ SimplyGo Mini comes with everything you need for an independent lifestyle. Spread the weight of the SimplyGo Mini for greater independence with the backpack accessory. Available separately in black or brown.




▲ Easy to use touchscreen control showing current setting and remaining battery life



◀ Battery options: The SimplyGo Mini can be fitted with either a standard or extended battery. Desktop charger also available separately.

Standard Battery Package

- SimplyGo Mini main unit
- Rechargeable standard battery
- Carry bag with shoulder straps
- Accessory bag
- AC mains power supply unit
- DC power supply unit
- User manual
- Quick start guide
- Nasal cannula



Approximate battery run-time
(in hours:minutes)

	Standard Battery	Extended Battery
Setting 1	4:30*	9:00*
Setting 2	4:30	9:00
Setting 3	2:30*	5:00*
Setting 4	2:00*	4:00*
Setting 5	1:30*	3:00*

*Battery run-times have been collated from various sources and are not officially published by the manufacturer. These run-times should be used as an approximate guide only.

*Run-time dependent on breathing rate

Extended Battery Package

As above including:

- Extended battery instead of standard battery

23.9 (h) x 21.1 (w) x 9.1 (d)
in cm with standard battery installed

Concentration	90% -3% /+ 6%
Weight	5.0lbs (2.3kg) - Std. batt. 6.0lbs (2.7kg) - Ext. batt.
Power	AC: 100-240v DC: 19v
Noise	42db(A) on setting 2 (pulse)
Batt. charge time	Standard: 4 hours Extended: 8 hours
Warranty	3 years. 1 year on sieve beds, 90 days on battery & accessories
Airline approved	Yes, subject to airline's terms & conditions
CPAP/BiPAP	This device may be suitable for use whilst sleeping. This device cannot be used with a CPAP or BiPAP device.

intermedical
HOMECARE

0800 028 2194

www.portableoxygen.co.uk

Intermedical (UK) Limited, Homecare Division
Unit 6 Mill Hall Business Estate, Aylesford, Kent, ME20 7JZ, United Kingdom

PHILIPS
RESPIRONICS
UK AUTHORISED
DISTRIBUTOR

globalgroup
ISO 13485
UKAS
MANAGEMENT
SYSTEMS
5965
Cert: 87M18032

Trusted in **Respiratory Health**