

To place your order please call
0800 028 2194

Or buy online at www.intermedicaldirect.co.uk/homecare

Intermedical (UK) Limited, Homecare Division
Unit 6 Mill Hall Business Estate, Aylesford, Kent, ME20 7JZ, United Kingdom

Trusted in Respiratory Health

PHILIPS
RESPIRONICS
SimplyGo



Optional accessories
for your Philips Respironics SimplyGo

intermedical
H O M E C A R E

Philips Respironics SimplyGo | Accessories List



Battery pack

Highly Recommended

Consider purchasing extra battery packs to extend your machines operation time. Up to 3 hours run-time on setting 2 pulse flow.

RE-1082662



Desktop charger

Highly Recommended

More than one battery? You may wish to consider a desktop charger. Charge one battery whilst using another.

RE-900-119



Liberty battery kit

Highly Recommended

Use two batteries at once to extend the overall battery operation time. Includes battery and external battery module.

RE-1116269



Carry bag and shoulder strap

Stylish protective carry bag available in brown

RE-1082663



Accessory bag

Keep your accessories all together. Perfect for when travelling.

RE-1083696



Humidifier kit

Humidifier bottle, connector tubing and pouch for SimplyGo

RE-SGHKIT



External battery module

Connect another battery to the SimplyGo using the external battery module

RE-1109643



Mains AC power supply unit

Replacement power supply unit that connects the machine to the mains power supply. Supplied with UK plug.

RE-1082661



DC power cable

Replacement power supply unit that connects the machine to a 12v socket

RE-1083692



DC cord for airlines

Power cable for use on selected aircraft. *Before purchasing, please contact your airline to ensure compatibility.*

RE-1083693



Cart with telescopic handle

Large wheel cart for easy transportation of the SimplyGo

RE-1074885